
CYTCO-K™

C106K – C107K – C108K – C109K – C110K – C219 – C117K – C0118K – C117K TK – C111P

**FOR DENTAL USE ONLY****GB****DIRECTIONS FOR USE CYTCO-K™****0) COMPOSITION**

Titanium

1) INDICATIONS FOR USE

These instruments / products are to be used only in a clinical or hospital environment, by qualified users.

Application for the product: In case of insufficient residual tooth substance (< 4 mm), the CYTCO-K is needed to support the coronal restoration.**2) CONTRA INDICATIONS**

Insufficient residual dentin: At least 2 mm of tooth structure is required around the post.

3) WARNINGS

To obtain improved distribution of forces on the dental reconstruction, natural dentine must be present over 2 or 3 mm all around the pivot

4) PRECAUTIONS

- The Post is a Single Use device.
- Re-use of Cytco-K can increase the risk of cross contamination, poor cementing and breakage
- The post must be sterilized before insertion in the canal.
- Avoid touching the posts with your fingers after cleaning.
- Dentsply Maillefer recommends the use of a rubber dam.

5) ADVERSE REACTIONS

In the present technical state, no adverse reaction has been reported so far.

6) STEP-BY-STEP INSTRUCTIONS

Root preparation and selection of root post

After filling the root canal, the appropriate CYTCO-K root post is selected by comparing the outlines on the transparent sheet (REF C 700K) with the radiograph of the tooth. The scale adopted for these drawings (1.2:1) corresponds approximately to the enlargement, obtained in dental radiography.

- Preliminary drilling: CYTCO-K Penetration Drill Ref C107K (Fig. 1)
- Precision Reaming and Counterboring: CYTCO-K Combination Drill Ref C108K (Fig. 2)
- Checking the reamed shaft and the future position of the Root Post: CYTCO-K Gauge Ref C109K (Fig. 3)
- Inserting and Cementing the CYTCO-K Root Post: using Wrench Ref C110K (Fig. 4)
- Restoring the tooth

References	Description	Packaging
C117K	Complete set	Single Unit
C117K TK	Trial kit	Single Unit
C106K	Posts	Blister pack cont. 12 pieces
C107K	Penetration drill	Box cont. 6 pieces
C108K	Combination drill	Box cont. 6 pieces
C109K	Gauge	Blister pack cont. 1 piece
C110K	Wrench	Blister pack cont. 1 piece
C111P	Sterilize Mount. block	Single unit
C0219	Mandrel	Blister pack cont. 1 piece
C700K	Transparent sheet	Single unit

7) DISINFECTION, CLEANING AND STERILIZATION

Reprocessing procedure

Foreword

For hygiene and sanitary safety purposes, all instruments not marked "sterile" must be cleaned, disinfected and sterilized before usage to prevent any contamination.

Area of application

Disinfection and sterilisation before usage concerning:

Implantable devices:

Dentinal and radicular posts made of steel, titanium and glass fibers.
Supports, kits and organiser systems for posts.

General recommendation

- 1 - Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use anticorrosion disinfecting and cleaning agents
- 2 - For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- 3 - The user is responsible for the sterility of the product on thee cycle as well as for the usage of damaged or dirty instruments where applicable after sterility.
- 4 - Single use marked instruments are not approved for re-use.
- 5- The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.

Step-by-step procedure

	Operation	Operating mode	Warning
1a.	Automated Cleaning with washer-disinfector	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments or posts. - Put them in the washer-disinfector (Ao value > 3000 or, at least 5 min at 90 °C). 	<ul style="list-style-type: none"> - Discard any posts with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the washer-disinfector use kits, supports or container. - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly.
OR			
1b	Manual Cleaning or assisted by an ultrasonic device	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments. - Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable. 	<ul style="list-style-type: none"> - No visible impurities should be observed on the posts. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Follow instructions and observe concentrations and time given by the manufacturer (see also general recommendations). - The disinfecting solution should be aldehyde free and without di- or triethanolamines as corrosion inhibitor.
2.	Rinsing	<ul style="list-style-type: none"> - Abundant rinsing (at least 1 min) 	<ul style="list-style-type: none"> - Use quality water in accordance with local regulations. - If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. - Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air.
3.	Inspection	<ul style="list-style-type: none"> - Inspect devices and sort out those with defects. - Assemble the devices (stops) 	<ul style="list-style-type: none"> - Dirty posts must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the posts.
4.	Packaging	<ul style="list-style-type: none"> - Place the devices in a kit, support or container to avoid any contact between instruments or posts and pack the devices in "Sterilisation pouches". 	<ul style="list-style-type: none"> - Avoid any contact between instruments or posts during sterilization. Use kits, supports or containers. - Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607.
5.	Sterilization	<ul style="list-style-type: none"> - Steam sterilisation at: 134 °C / 273°F-during 18 min. 	<ul style="list-style-type: none"> - The posts must be sterilized according to the packaging labelling. - Use only autoclaves that are matching the requirements of EN 13060, EN 285. - Use a validated sterilisation procedure according ISO 17665 - Respect the maintenance procedure of the autoclave device given by the manufacturer. - Use only this recommended sterilization procedure. - Control the efficiency (packaging integrity, no humidity, colour change of sterilisation indicators, physico-chemical integrators, digital records of cycles parameters). - Traceability of procedure records
6.	Storage	<ul style="list-style-type: none"> - Keep devices in sterilization packaging in a dry and clean environment 	<ul style="list-style-type: none"> - Sterility cannot be guaranteed if packaging is open, damaged or wet. - Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period).

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