



# RADIX-ANKER® POSTS

C0193, C193T / C195C, C195CT, C195D, C195DT, C195ET  
C0261, C0310, C310CT, C0320, C0322

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## FOR DENTAL USE ONLY

### DIRECTIONS FOR USE

0) COMPOSITION: Stainless Steel or Titanium.

### 1) INDICATIONS FOR USE

These products have to be used only in hospital environments, clinics or dental offices by qualified dental personnel.

**Application for the product:** In case of insufficient residual tooth substance (< 4 mm), the Radix-Anchor® is needed to support the coronal restoration.

### 2) CONTRA-INDICATIONS

Insufficient residual dentin: at least 2 mm of tooth structure is required around the post.

### 3) WARNINGS

To obtain improved distribution of forces on the dental reconstruction, natural dentine must be present over 2 or 3 mm all around the pivot

### 4) PRECAUTIONS

The post must be cleaned with alcohol before insertion in the canal. Avoid touching the posts with your fingers after cleaning. Dentsply Mallefer recommends the use of a rubber dam.

### 5) ADVERSE REACTIONS

In the present technical state, no adverse reaction has been reported so far.

### 6) STEP-BY-STEP INSTRUCTIONS

The following rules must be respected in preparing the anchorage cavity and in seating and inserting the Radix-Anchor:

1. The crown part of the canal will remain open or will be opened by manual instruments down to penetration depth of the post

2. The canal cavity is initially drilled with the non-cutting tip penetration drill C0220.

3. The seating for the Radix-Anchor head is drilled with the root facer C0217 (The stability of Radix-Anchor® is thus ensured).

4. Precision drilling is performed manually with the spiral bur C0216 and the mandrel C0219.

5. Using the gauge C0191 the preparation for the Radix-Anchor® can be checked as well as its position. The cutting of threads in the root cavity will be realized by screwing the post in place alternately by a half-turn clockwise and a quarter-turn anti-clockwise.

6. Unscrew the post to insert the sealing cement in the root canal and on the post itself.

7. When the anchor post is reinserted, it must first be turned to the left until it engages with the pre-cut thread, and only then screwed to the right.

8. Clean the head and root section, then apply a dentine adhesive. The composite material, resin or amalgam is applied for crown restorations, temporary crowns and fillings.

### Safety and Liability

Dentsply Mallefer declines any liability or compensation for possible damage due to:

- use of instruments which do not belong to the system and which may impair the function.
  - failure to observe the instructions for use.
- The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the directions for use.

### 7) DISINFECTION, CLEANING AND STERILIZATION:

See general instructions in appendix

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## GB REPROCESSING PROCEDURE FOR DENTAL INSTRUMENTS AND IMPLANTABLE RADICULAR DEVICES

### Foreword

For hygiene and sanitary safety purposes, all instruments not marked "sterile" must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent ones.

### Area of application

Disinfection and sterilisation before first usage and reprocessing procedures concerning:

#### A1 Instruments:

Cutting instruments, (hand and engine driven) such as:

- Endodontic instruments (files, broaches, reamers, enlargers, endodontic burs, ultrasonic inserts);
  - Rotary cutting instruments (Diamond burs, tungsten carbide burs, stainless steel drills, carbon steel burs).
- Root canal filling instruments (Pluggers, spreaders, compactors);  
Supports, kits and instrument organisers  
Hand instruments and clamps.

#### A2 Implantable devices:

Dentinal and radicular posts made of steel, titanium and glass fibers.  
Supports, kits and organiser systems for posts.

#### A3. Contra angle

B. Filling material: Only chemical disinfection (no sterilisation)

General recommendation

- Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use antirust-corrosion-disinfecting and cleaning agents.
• For your own safety, please wear personal protective equipment (gloves, glasses, mask).
• The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.
• Limitations and restrictions on reprocessing:
- The individual DFU indicates if the useful life of a device might be reduced by the number of reprocessing cycles. Furthermore, the appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of colour coding or marking, are indications that the devices are notable to fulfil the intended use with the required safety level.
• Single use marked instruments are not approved for re-use.
• The water quality has to be convenient to the local regulations especially for the last rinsing step with a washer-disinfector.
• Tungsten carbide burs, plastic supports, hand instruments and NiTi instruments are degraded by Hydrogen Peroxide (H2O2) solution.
• NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5%.
• Instruments made of aluminium are degraded in presence of caustic soda solutions with mercury salt. Do not use acid (pH < 6) or alkaline (pH > 8) solutions.
• The washer-disinfector is not recommended for instruments made of aluminium, tungsten carbide or carbon steel.

Step-by-step procedure

A. Devices

Table with columns: Operation, Operating mode, Warning, and a grid of checkboxes for various procedures like Disassembly, Pre-Disinfection, Rinsing, Automated Cleaning, etc.

B. Filling material

Table with columns: Operation, Operating mode, Warning, and a grid of checkboxes for Disinfection.

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